

LANDAUER OPTIMIZE™



Data in our hands. Insight in yours.
A physicist at your fingertips.



Your dose management solution should provide

- ✓ Industry leading expertise
- ✓ Daily monitoring
- ✓ Best practice guidance
- ✓ Optimization for quality care

OPTIMIZE: On-demand physics expertise via a simple subscription



Insight

Alert Review & Analysis

Let us take on the burden of managing alerts so issues are identified and corrected immediately



Expertise

Best Practice Guidance

Recommendations for balancing dose and image quality, policy and procedure development and standardizing care enterprise-wide



Compliance

Benchmarking & Protocols

Hands-on protocol management, size-specific and pediatric protocols, CT and fluoro benchmarking and recommendations for dose index ranges



Simplicity

Automated Data Collection


Daily insights from data captured using our cloud-based software

Optimization at Work

Daily monitoring and unlimited access feed a cycle of improvement helping you to stay compliant and deliver quality patient care

Daily Monitoring

Your dedicated LANDAUER physicist performs their routine review of your alerts and performance



Insights in Your Inbox

Receive an email with their findings and recommendation



Expertise




Opportunities to optimize are identified



Unlimited Access

Click to review the insight in the system or to schedule a meeting to discuss




Confidence in your daily practice

The longer a protocol goes without being optimized or incidents go unaddressed, the greater risk to patients. OPTIMIZE is the answer.


Quarterly Reviews

Without daily governance issues can go unidentified



Daily Monitoring

Faster resolutions with OPTIMIZE by our physics experts we reduce outliers quickly for improved patient care



Time-Savings

To alleviate the burden on your team, we analyze your data and share only what you need to act on




Getting Connected

Here's how it's done

- Install an on-premise site server (VM)
- VM connects your PACS to our cloud-based software
- Real time analysis begins

Here's how we make it easy

- No PHI leaves your hospital network
- Eliminates need to export data
- Simplifies upgrades



Why LANDAUER?



5 Million

Exams
analyzed
each year

2000
Insights annually

100% Compliance

250 Customers
Leveraging our expertise



20,000

Protocols reviewed
each year

Visit landauer.com/OPTIMIZE

LANDAUER®
OPTIMIZE