


# Enterprise Dose Dashboard

Available via myLDR.com

## Gain access to enterprise-wide dose dashboards to compare multiple facilities across an entire organization

Enterprise Dose Dashboards enable you to more effectively track dose levels across an enterprise. At the same time, you can more easily track improvements as a result of conducting training, investing in shielding, or buying other protective equipment. The dashboards provide data to help drive change and improve radiation safety across the enterprise. Also, see how your facilities stack up against each other and the LANDAUER ALARA Index\*. Dashboards are provided for ALARA Account Comparison, ALARA Subaccount Comparison and Pregnant Participants Review (Cumulative and Monthly). All dashboards are also accompanied by data in a sortable table format.

### ALARA Account Comparison - Across the enterprise by Account



Track and identify ALARA 1 and 2 performance and compare it across Accounts. Use filters to select from the last eight quarters of data.


Utilize the dashboard to see how Accounts compare to the ALARA 1 and 2 enterprise average and Landauer ALARA Index\*.

\*LANDAUER ALARA Index represents all clients' participants who have exceeded the levels. Calculation is based upon the previous three calendar years of data.


### ALARA Subaccount Comparison - Across the enterprise by Subaccount

With a more targeted view, easily identify and compare selected Subaccounts across the enterprise with ALARA 1 and 2 exceeders.

This is especially helpful for comparing similar Subaccounts at different Accounts where you know or suspect the dose is typically higher. For example, filter for all Cath Labs across the enterprise to see where there are opportunities for training or improvement at the Subaccount level.



## Pregnant Participants Review - Cumulative View



Use the first dashboard to visualize the current number of active declared pregnant workers across Accounts.

## NEARSHORE HEALTHCARE Cumulative Fetal Dose vs Pregnancy Duration

Data As Of  
03-19-2015



Each dot on the dashboard above represents a pregnant participant.

Use the second dashboard to compare Cumulative Fetal Dose over Pregnancy Duration by participant. It provides insights for dose and duration in the following ways:

- where participants are trending towards the 500 mrem limit for the pregnancy
- where participants are over the limit
- where a participant is still being identified as pregnant (beyond the expected 9 month duration)

## Pregnant Participants Review - Monthly View



Visualize the number of active declared pregnant workers broken down by dose ranges across Accounts.

There are reference lines at 40 mrem and 50 mrem to quickly identify participants with fetal doses at or greater than those values.

Use filters to select from the last 12 months of data.

Contact your Regional Sales Manager or email [myLDRsupport@landauer.com](mailto:myLDRsupport@landauer.com) to learn more about how your facility can benefit from Enterprise Dose Dashboards.